

BABS1202 PRACTICAL EXAM RUBRIC

Tasks	0 marks	1 mark	2 marks
Task 1 Gram-staining	Uneven distribution of bacterial cells AND incorrect Gram status	Uneven distribution of bacterial cells OR Incorrect Gram-status for stained sample.	Even distribution of bacterial cells AND Correct Gram-status for stained sample.
Task 2 Light microscopy	Student could not independently focus on a field of view AND incorrectly identified Gram-status and shape of visualized sample.	Student could not independently focus on a field of view OR incorrectly identified Gram-status and shape of visualized sample.	Student independently focused on a field of view AND correctly identified Gram-status and shape of visualised sample.
Task 3 Streak-plating	Absence of any dilution in microbial growth via streaking.	Successful isolation of single colonies following streaking from the primary inoculum.	-
Task 4 Sterile broth transfer	Presence of turbidity (contamination) in BOTH bottles following incubation. Absence of turbidity in inoculated bottle following incubation.	Presence of turbidity in EITHER bottle following incubation.	BOTH bottles remain sterile following broth transfer and incubation.
Task 5 Broth culture inoculation	Absence of turbidity in inoculated bottle following incubation.	Presence of turbidity in inoculated bottle following incubation.	-
Task 6 Serial dilution: agar plate (viable plate count)	Number of colonies on final dilution plate differs by more than an order of magnitude compared to predicted number.	Number of colonies on final dilution plate within same order of magnitude as predicted number.	-
Task 7 Serial dilution: 96-well plate	Incorrect 1:10 serial dilution of food dye across 8 wells.	Correct 1:10 serial dilution of food dye across 8 wells.	-

Total marks available = 10